

Remote Operated Obstetric OT Table

Model: LM-OTGP



Lumens Medical Pvt. Ltd.
Life Enhancing Technology



Optional:
Lateral Tilt

FEATURES:

- Electromatic drives with noiseless, jerk less and smooth positioning.
- Three section upholstered table top ensures maximum comfort to the patient.
- Table positions i.e. Height adjustment, Trendelenburg, Reverse Trendelenburg, Back Section/Chair Position are controlled by corded handset with graphical icons.
- Detachable leg sections offer perineal access & allow to attach drainage tray.
- Integrated weight compensation by manual adjustment of leg section.
- Large size stainless steel drainage tray
- Non-corrosive stainless steel 304G construction including table top frames, accessories rail, base cover and column telescopic covers.
- Easy for cleaning & hygiene.
- The powered movements are activated by Imported Linear Actuators.
- Four easy glide, antistatic castors, two of which are swivelling, ensure mobility.
- Pedal operated central braking system secures the table in its place while in use.

TECHNICAL SPECIFICATIONS:

- Top dimension: L 1750 x W 600 mm
- Height adjustment with pad: 700 mm – 950 mm
- Trendelenburg / Reverse Trend: 15° / 15°
- Back Rest (upward / downward): 60° / 0°
- Leg Rest (upward / downward): 15° / 90°
- Power supply: 24VDC
- Patient Weight Capacity: 150 kgs (330 lbs)



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

An ISO 9001, ISO-13485, CE, CDSCO CLASS A&B Certified Company
2-24-113/3/B/NR, Laxmi Narayana Nagar, Uppal, Hyderabad.

Tel.: 9849466520, 7997192056

Email : lumensmedical@gmail.com www.lumensmedical.com



Height adjustment



Trendelenburg / Rev. Trend.



Back Section adjustment



Lithotomy Position