

Germans call it
Blitz
We call it
Pure Natural
Surgical Light



Lumens Medical Pvt. Ltd.™
Life Enhancing Technology

Blitz Series Hi Depth Surgical Illumination Systems

POWERED BY HIGH DEPTH LED'S

"Blitz L6 Single Dome" is a High depth Surgical illumination System with total shadow dilution. The low glare illumination with high depth enables Surgeon to do complicated & time consuming Surgeries without any Eye strain.



Model: **Blitz L6 Single Dome**

Key features:

- Focus adjustment available
- Colour Temperature 3500 - 5000 K
- Intensity Adjustment 10% - 100%
- CRI >95
- Display clear LCD Screen
- LED Life 60,000 Hours



Intensity Controllers

Remote sensor

White Lights

Yellow Lights

Removable Sterilizable Handle



SPECIFICATIONS

- Light Intensity: 1,40,000 Lux
- Endo Mode available
- Colour Rendering Index (CRI): 95
- Colour Temperature 3500 - 5000 K
- Life span of LED 60000 Hours
- Round shaped lamp head.
- Dome: ABS material
- Control panel: manual
- Laminar flow compatible.
- Illuminated field diameter - 180mm - 300 mm adjustable
- Illumination Depth - 100 to 140 cm.
- Dome Diameter - 500 mm
- Variation with Digital display 10-100% in 10 steps
- User selectable intensity
- On / Off Switch available
- Scratch proof glass with interchangeable lens
- Spring arm with 360° rotation
- User friendly Sterilizable handle at Center of Lamp Head
- Power supply 220-240V - 50 Hz Single phase
- CDSCO Certified, ISO 13465, ISO 9001, IEC 60601-1, CE

Our Pride

Over 2000 surgeries
Performed under our
Surgical lights per day

Installation base

More than 1000 satisfied
surgeons using Lumens
Surgical lighting
systems all over India

Quality Certifications

ISO 9001
ISO 13485
IEC 60606-1-1
IEC 6060-1-1-2
CDSCO, CE



**MADE IN INDIA
FOR THE WORLD**



Lumens Medical Pvt. Ltd.™

Life Enhancing Technology

An ISO Certified, ISO-13485, CE, IEC 60601-1, IEC 60601-1-2, CDSCO CLASS A&B
2-24-113/3/B/NR, Laxmi Narayanna Nagar, Uppal, Hyderabad.

Tel.: 9849466520, 7997192056

Email : lumensmedical@gmail.com www.lumensmedical.com